Country Code: _____

Student Code: _____

The 21st INTERNATIONAL BIOLOGY OLYMPIAD

11th – 18th July, 2010

Changwon, KOREA

PRACTICAL TEST 3

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GENETICS AND CELL BIOLOGY

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Total Points: 50

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Duration: 90 minutes

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Dear Participants, 敬愛的參賽者

- ☺ In this test, you have been given the following 2 tasks: 本試卷包含兩大題

Task I (35 points)

- (1) Study of promoter-driven regulation of gene expression.** (20 points)

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- (2) Characterization of the relationship between genotypes and phenotypes**

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Task II: Observation of meiotic cells in preserved rye anthers (15 points)

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- ☺ Write down your results and answers in the **Answer Sheet**. **Answers written in the**

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- ☺ Please make sure that you have received all the materials listed for each task. If any of

the listed items is missing, please raise your hand. 確認所有材料的種類與數量都與試題

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- ☺ Stop answering and put down your pencil **immediately** after the end bell rings. The

supervisor will collect the Question Paper and the Answer Sheet. 鈴聲響時，立刻停止

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Good Luck!! 助你好運

GENETICS AND CELL BIOLOGY

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This practical test is composed of 2 tasks.

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TASK I. (35 points)

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(1) Study of the promoter-driven regulation of gene expression

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(2) Characterization of the relationship between genotypes and phenotypes

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On individual Table

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1. Fluoro-spectrophotometer

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2. Microfuge tubes containing 50 μ L each of nine differently-labeled plant extracts. Two identically labeled tubes are provided for each type of extract ($2 \times 9 = 18$ tubes) The transparent tubes are for the protein assay, and the black tubes are for fluorescence measurements.

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3. 12 mL Bradford reagent in a 15 mL plastic tube (Bradford reagent is used to determine concentration of protein)

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4. 1 mL of 1 mM MUG (fluorescence substrate to measure GUS activity) in a microfuge tube

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5. 12 mL of stop reagent for the GUS (enzyme β -glucuronidase which converts MUG into MU) reaction in a 15 mL plastic tube

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6. Two DNA size-marker tubes (labeled M, 50 μ L each) and eight tubes containing *Eco*RI-digested DNA (labeled P1~P8, 50 μ L each)

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7. Two microfuge tubes labeled as GUS BL and Pro BL, respectively.

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8. Three micropipettes (one each for 10-100 μ L and 100-1000 μ L, and a fixed volume pipette for 20 μ L)

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9. A box of yellow tips for the 20 μL and the 10-100 μL micropipettes

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running buffer. If your gel is broken, raise your hand for assistance.

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12. A tip disposal container 微量吸管放置架。

13. Polygloves 橡膠手套。

14. 25 cuvettes for the Fluoro-spectrophotometer

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16. A timer 計時器

17. A Scotch tape 膠帶

18. An ice bucket filled with ice 裝有冰塊的冰桶

19. Microfuge tube racks 微量離心管管架

20. Green card 綠色卡片

On the common equipment table

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
1. Gel documentation system equipped with a UV source

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Handling of Micropipettes

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Adjustment method 調整模式

Turn the plunger to set the volume to the desired value, which can be seen in the display window. Remember that each micropipette has designated range of volumes as indicated on the pipette. Do not exceed the limits of this range.

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Usage method 使用方法

- 1) Secure the pipette tip to the tip holder. Gently push down the plunger to the first stop.

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- 2) Hold and lower the tip down into the solution to a depth of 2~4 mm. Release the plunger slowly to allow it to return to its original position.

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- 3) Remove the pipette from the liquid, and transfer the contents to the desired tube. Push the plunger to the first stop and then push further to the second stop to discharge the solution completely from the tip.

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
- 4) Remove the pipette from the tube and release the plunger. Eject the used tip into the tip disposal container by pressing the tip-ejector.

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Operating Instruction for the Fluoro-Spectrophotometer (measures fluorescence of MU and absorbance of proteins at 595 nm)

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A: Cuvette holder for protein measurement

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B: Cuvette holder for fluorescence of MU measurement

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FCT: Function key

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BL: Blank key

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TS: Test sample key

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PWR: Power key

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Usage method 使用方法

Important: Please be sure not to touch the light path of cuvettes.

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- 1) Press the PWR () button to turn on the machine. The display window will be turned on after a beep.

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- 2) To set the blank sample to zero, insert the blank cuvette in an appropriate holder (use cuvette holder A to measure protein concentration, and cuvette holder B to measure GUS activity). The cuvette indicator will be turned on ( for the holder A and  for the holder B).

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Note: Two blank samples for measurement of GUS activity and amounts of proteins are provided in the microfuge tubes labeled as GUS BL and Pro BL, respectively.

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- 3) Press the BL button, and the blank indicator () will appear when the blank is set at 0.0.

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- 4) To measure a sample, remove the blank cuvette and insert the test cuvette in the same cuvette holder, and press the TS button. The result will be displayed after 5-10 seconds, and the

indicator will appear in the display window ()

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- 5) To end the machine, keep the PWR button pressed till beep is heard.

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Operating Instruction for the DNA Gel Electrophoretic Apparatus

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- 2) After verifying that the operation switch of the power supply is OFF, close the migration tank lid.

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Do this as follows:

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
- (1) First, insert the 2 tabs on the cover into the holes in the migration tank.

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
- (2) Then, rotate the cover forward to close it.

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
- 3) Set the voltage to “Half” using the output selection switch.

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- 4) Push the operation switch to start the migration.

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- 6) In this experiment, the gel running time should be 30 min. Make sure to turn the operation switch OFF when the running is finished

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Using the gene X-fused GUS reporter gene to analyze hormonal effects on gene expression and to characterize the hormone-responsive elements in the promoter.

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Plants respond to their hormones by regulating hormone-responsive genes. Within a gene promoter, a specific DNA sequence(s), the *cis*-element, dictates the proper time and amount of gene expression. Regulation is primarily controlled by a hormone-responsive transcription factor(s) that binds specifically to this region, resulting either in gene activation or suppression.

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
In this task, you will examine the mode of hormonal regulation in the hormone-responsive gene X of *Arabidopsis*. To find the hormone-responsive regions, in the promoter and to understand the mode of hormonal regulation of gene X expression, the promoter of gene X is divided into A~C (each of these domain may function as enhancer, silencer or minimal promoter).

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Then, a variety of *Arabidopsis* transgenic plants expressing the GUS (β -glucuronidase) reporter gene under the control of the different regions of the promoter, as diagramed below, was generated. The GUS will be produced when the promoter of gene *X* is activated. The GUS enzyme converts MUG into MU, and its activity can be measured by quantifying MU fluorescence using a fluoro-spectrophotometer.

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Q1. The purpose of the first experiment is two-fold: (1) to find the promoter region containing a hormone-responsive *cis*-element and (2) to investigate the effects of different hormone H concentrations on gene X expression. All transgenic plants (WT, dA, dAB, and dABC) were treated with either 1 μ M or 100 μ M of hormone H. To assess the level of GUS expression, plant extracts were prepared from these treated plants. (See the table in the materials and method section.)

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Using the methods described in the next section, measure the fluorescence value and absorbance at 595 nm of each 50 μ L plant extracts. Based on these measurements, calculate the amount of MU (nmole MU/50 μ L plant extracts), the amount of proteins (μ g/50 μ L plant extracts), and the resulting GUS activity (nmole MU/ μ g protein/min) for each extract. Record your results in Table 1 in the answer sheet to find answers for **Q1.1**, **Q1.2**, and **Q1.3**.

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Measurement of fluorescence and determination of MU amount

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- 1)-1. Turn on and set the fluoro-spectrophotometer to zero with 500 μ L of the blank sample labeled GUS BL.

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- 1)-2. Take a microfuge tube of plant extracts (each tube contains 50 μ L extracts) prepared from each WT-O or hormone-treated transgenic plant, and **mix well** (by gentle tapping) with 50 μ L of 1 mM MUG solution. Start with the labeled WT-O and proceed in an order shown in the table in Materials and Equipment.

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- 1)-4. Stop the reaction by adding 900 μ L of stop reagent (1M sodium carbonate in GUS extraction buffer) into each 100 μ L reaction solution In the same order you added MUG. **Mix well** by tapping.

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- 1)-5. Take 500 μ L of the finished mixture from each tube, and measure the fluorescence using the fluoro-spectrophotometer.

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- 1)-6. Calculate the amount of MU in the sample using the formula provided below. Record the fluorescence value and the calculated amount of MU in Table 1 in the answer sheet. This is the amount of MU produced from each of 50 μ L plant extracts.

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$$Y = 0.04 X + 2.5$$

Y: the amount of MU (nmoles 50 μ L plant extracts)

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X: the measured fluorescence value [from step 1)-5]

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Measurement of absorbance at 595 nm and determination of protein amount

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- 2)-1. Turn on and set the fluoro-spectrophotometer to zero with 500 μ L of the blank sample labeled Pro BL.

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- 2)-2. Take a microfuge tube with extracts (each tube contains 50 μ L extracts) prepared from each WT-O or hormone-treated transgenic plant, and mix well with 950 μ L of Bradford reagent. Incubate at room temperature for 5 min.

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2)-3. Take 500 μL of the reaction mixture from each tube, and measure the absorbance at 595 nm using the fluoro-spectrophotometer.

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2)-4. Calculate the amount of proteins using the formula provided below. Record the absorbance at 595 nm and the calculated amount of proteins in Table 1 in the answer sheet. This is the amount of proteins contained in each of 50 μL plant extracts.

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$$Y = 98X + 2.8$$

Y: the amount of protein ($\mu\text{g}/50 \mu\text{L}$ plant extract)

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X: the measured absorbance at 595 nm of the solution [from step 2)-3]

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Calculation of GUS activity GUS 活性計算

3)-1. Considering that this GUS enzyme reaction was performed for 10 min [refer to 1)-3], calculate GUS activity in nmole MU/ μg protein/min and record the value in Table 1 in the answer sheet.

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Table 1 is worth of 9 points.

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Q1.1. (4 points)

Based on your results in <Table 1>, put a checkmark (✓) in the appropriate box of each plant in Table **Q1.1** in the answer sheet.

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Note: - stimulation: more than **3-fold** increase in gene *X* expression

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- no effect: less than **3-fold** increase in gene *X* expression

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Q1.2. (6 points = 2 × 3)

Based on your previous conclusions in **Q1.1**, determine the regulatory function (enhancer, silencer, or minimal promoter) of each *cis*-element (A~C). Put a checkmark (✓) in the appropriate box in Table **Q1.2** in the answer sheet.

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Q1.3. (1 points)

How does 100 μM of hormone H regulate the expression of gene X? Based on your finding from <Table 1>, determine the mode of action of hormone H. Put a checkmark (\checkmark) in the appropriate box in Table **Q1.3** in the answer sheet.

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
Part II. (15 points)

A co-relationship analysis between genotype and phenotype, and the prediction of gene pool frequencies using Hardy-Weinberg mathematics.

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Q2. Gene *Y* encodes a protein that regulates plant growth. The schematic figure below depicts the region of gene *Y* in genomic DNA and a point mutation within.

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There are eight plants with homozygous (YY or yy) or heterozygous (Yy) genotype, showing either wild type or dwarf phenotypes (Y : wild type allele, y : mutant allele. The alleles Y and y do not specify whether they are dominant or recessive). To analyze the genotype of these plants, the 1 kb region of gene Y was amplified by PCR. This fragment was then digested with *EcoRI* restriction enzyme, which cuts GAATTC sequence. Other than the *EcoRI* site created by the point mutation, there is no other *EcoRI* recognition sequence in gene Y . Using the protocol described below, perform a gel electrophoresis of the *EcoRI*-digested PCR products.

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Genotyping of gene Y by gel electrophoresis

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Note: Always wear polygloves during the experiment !!!

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- (1) A total of ten microfuge tubes are provided: two DNA size marker tubes (M) and eight tubes containing *Eco*RI-treated PCR product from Plants 1~8 (P1~P8, respectively).

Starting from left in the order of M, P1 – P8, M, load 20 μ L out of 50 μ L DNA solution into each well of a prepared agarose gel in the electrophoresis apparatus. Use the 20 μ L micropipette to load samples. Change pipette tip for each sample.

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Note: The DNA size marker solution contains 0.4, 0.6, and 1.0 kb DNA fragments.

DNA loading buffer and DNA-staining dye are already included in each tube.

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- (2) Refer to <Operating instructions for DNA gel electrophoretic apparatus> to put the cover on the electrophoresis apparatus, to turn on the apparatus, and to run the electrophoresis.

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Note: Upon starting the electrophoresis, make sure that the output indicator LED is lit and that bubbles are forming on the platinum electrodes.

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*** IMPORTANT: While the gel is running, proceed to TASK II !!!**

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- (4) Turn off the apparatus. Then, raise the **green card** to request help for photography of the agarose gel.

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Note: The assistant will bring a gel transfer box to you. Make sure that your student code is on the box.

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- (5) When you receive the agarose gel picture, attach it to Q2.1 of the answer sheet using Scotch tape. Label the number of each plant (P1~P8) on each lane of the gel picture.

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- (6) In Table Q2.2 in the answer sheet, put checkmarks (✓) to designate the size of DNA fragments and the genotype of each plant.

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Q2.1. (3 points)

Attach the agarose gel picture to a space given on the answer sheet. And label the number of each plant (P1~P8) on each lane of the gel picture.

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Q2.2. (4 points)

Determine the size of DNA fragment(s) and the genotype (*YY*, *Yy* or *yy*) of each plant. Put a checkmark (✓) in the appropriate box in Table **Q2.2** in the answer sheet.

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Q2.3. (2 points)

Based on the genotype and phenotype of each plant given in **Q2.2**, deduce the characteristic of the mutation. Put a checkmark (✓) in the appropriate box in the Table **Q2.3** in the answer sheet.

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Q2.4. (2 points)

If you cross Plant 1 with Plant 3 (from **Q2.2**), what is the probability (%) that an offspring will be a dwarf plant? Write your answer in the answer sheet.

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Q2.5. (4 points)

The eight plants in Q2.2 represent a population. If this population produces 10,000 plants in the next generation, what would be the expected number of heterozygous and dwarf offspring, respectively? (Assume that this population is in Hardy-Weinberg equilibrium.).

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TASK II. (15 points)

Observation of meiotic cells in preserved rye anthers

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Background 背景說明

Using a light microscope, you will observe meiotic cells in preserved rye anthers. Anthers at a specific stage of meiosis were selected and were preserved in 70% ethanol.

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Requirements – Overview

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Using the microscope, identify anther cells undergoing meiosis. In the space given in the answer sheet, sketch an image of meiotic cell you observe at 400X magnification (**Q3.2**)

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Procedure 步驟

- 1) Before you start observation, check for the presence of two small preserved anthers in the vial.

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
- 2) Take out the ceramic tile out of the tray, and put one glass slide on it.

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- 3) Observe your specimen under the microscope at 100X magnification, and find at least one cell undergoing meiosis. Then, observe **one** cell at 400X magnification and draw this image in the given area of the answer sheet (Q3.2). Make sure that this cell is at the center of your field of view. After you finish the drawing, raise the **red card**. The lab assistant will come to you and will take a photograph of the slide.

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1. In step (1), if the anthers won't come out, put the solution back into the vial using the disposable plastic pipet and repeat step (1) .

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2. Be careful not to break the anther in step (2)

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3. You may use a filter paper to remove excess 70% ethanol in step (3).

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4. Do not press too hard, or you may break the cells and/or the cover slip in step (7).

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5. You are provided with two anthers to prepare your specimen. If you fail to make good specimen with the first anther, please repeat the procedure and make another preparation using the other. However keep in mind that the time for your experiment is limited.

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
Q3. Answer the following questions.

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***Important:** You will see two types of cells under the microscope as shown in Figure Q3. The circled ones are examples of cells undergoing meiosis, and the rest are cells of the anther wall.

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400X

Figure Q3. Examples of cells undergoing meiotic cell division observed under a microscope.

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Q3.1. (1 point)

What kind of cells in the anther undergoes meiosis? Put a checkmark (✓) in the appropriate box in the answer sheet.

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Q3.2. (8 points)

Draw **one cell** undergoing meiosis at 400X magnification in the answer sheet. Do not label the drawing.

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Important : This cell must be at the center of your field of view when the picture is taken.

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Q3.3. (4 points)

At what meiotic stage are the cells? Put a checkmark (✓) in the appropriate box in the answer sheet.

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Q3.4. (2 points)

What is the amount of DNA in the cell undergoing meiosis that you observed and a cell of the anther wall, respectively? Put checkmarks (✓) in the appropriate boxes in the answer sheet.

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