PRACTICAL EXAM 2 MOLECULAR BIOLOGY

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GENERAL INFORMATION

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Max. total points 100

Exam duration 90 minutes

15 questions

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The exam consists of three parts.

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Part 1. Restriction enzyme map of plasmid (48 points)

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Part 2. PCR-based genotyping of yeast mutants (37 points)

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Part 3. Amino acid auxotrophy of mutant yeast (15 points)

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We suggest you read the entire exam file before you begin the lab work. In order to accomplish the entire exam, you have to run the three parts in parallel (see Figure below).

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
Figure 0.1: Suggested time table for completion of the three parts of the exam.

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1 Check materials and equipment

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Task 1

- 2 Setup restriction enzyme digestion
- 3 Incubation of digestion reaction
- 4 Setup agarose gel
- 5 Prepare samples and load on gel
- 6 Gel electrophoresis*
- 7 Take a picture of the gel + Questions

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Task 2

- 8 Setup PCR reactions
- 9 PCR reaction (in thermal cycler)
- 10 Prepare samples and load on gel
- 11 Gel electrophoresis*
- 12 Take a picture of the gel + Questions

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Task 3

- 13 Observe and deduce + Questions

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**Use the same gel, but different wells for the samples*

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INTRODUCTION

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MATERIALS AND EQUIPMENT

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In order to do your lab work, you need the materials A-V listed below. Please, ensure that these materials are available to you. **If anything is missing, contact the exam personnel by raising your pink card immediately – and no later than 15 minutes after the beginning of the exam.**

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Please notice

1. All liquids are provided in 2X excess amount of what is needed in order to do the analyses. Additional materials (including tubes and agarose gel) cannot be provided in case of spillage or errors during the set-up of experiments.
2. Remember, in order to collect liquids at the bottom of a tube, use the technique that you were shown on Monday.
3. Remember, one of the purposes of the exam is to test your lab work skills.
4. Tubes labeled with a blue line on the label should be **kept on ice at all times**.
5. All solutions are frozen so you will have to **thaw and mix them** before use.

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D	1	box of pipette tips 20-200 μ L 2-200 μ l 滴管尖	1, 2
E	1	rack for PCR tubes PCR管架	1, 2
F	1	rack for 1.5 ml tubes 1.5 ml 管架	1, 2
G	5	1.5 ml tubes (in envelope) 1.5 ml 小管(大信封中)	2
H	2	strips of 0.2 ml PCR tubes (in envelope). Do not break the tubes apart 0.2 ml PCR小管聯條(大信封中)，勿將這些管聯條分開。	1, 2
I	1	White bag with agarose gel with 12 wells. (Do not open the bag until you are ready to load the gel) 白色塑膠袋內有12孔電泳膠(直到你要注膠時，才打開袋子)	1, 2
J	1	OneRun electrophoresis system with running buffer 電泳系統及電泳緩衝液	1, 2
K	1	marker pen 麥克筆	1, 2
L	1	stopwatch/timer 計時器	1, 2
M	1	paper tissue 紙	1, 2
N	1	set of plastic gloves (handed out in the waiting room) 塑膠手套(在等待室發給)	1, 2
O	2	zip-lock bags (in envelope) 夾鏈袋(大信封中)	1
		Booklet with photos of yeast	

P	1	grown on eight different media 酵母菌在8種培養基的生長照片冊	3
Q	1	Chart of biochemical reactions for formation of various amino acids (on the wall of the workspace) 實驗台牆上的胺基酸生化合成路徑圖	3
R	1	1 Kb Plus DNA ladder (in envelope) 1 kb plus 表尺DNA照片((大信封中))	1, 2
S	1	Pink card for contact with exam personnel 通知試務人員的粉紅卡片	1, 2, 3
T	2	Small ID tags 小ID標籤	1,2
U	1	Big ID tags 大ID標籤	1
V	1	Water 540 μ L 540 μ L水	1,2

Figure 0.2: Material H: Tubes and ID label filled out with country code and student number.

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Ice bucket containing 1.5 ml tubes with a set of liquids.

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1. RESTRICTION ENZYME MAP OF PLASMID (48 POINTS)

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An experiment was designed to determine the subcellular localization of protein Alb in *Saccharomyces cerevisiae*. The experimental strategy required that the gene *alb* (999bp) was fused with the sequences encoding two different fluorescent proteins. This was achieved by cloning the gene into two plasmid backbones (pX and pZ) resulting in plasmid pX:alb and pZ:alb.

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E. coli was transformed with the two ligation mixes. Plasmids from two resulting *E. coli* transformants were purified, yielding DNA found in the tubes "Plasmid tube 1" and "Plasmid tube 2". Your task is now to determine which two of the four possible plasmids (pX, pX:alb, pZ or pZ:alb) are contained in each of the two tubes (one in each tube)

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
Node Id: b494ba5c17b6eddea97b2778

Indicate here by a tag, the presence of expected sizes of DNA fragments obtained by a complete PstI digestion of the four plasmids

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Table 1.1. Restriction enzymes and their optimal reaction conditions.
Table 1.1 各種限制酶的最適反應條件

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RESTRICTION ENZYME DIGESTION OF THE TWO PLASMIDS

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Q. 2

CHOICE OF RESTRICTION ENZYME (4 POINTS)

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Node Id: 2a3e13184a7a07ed1aff6caf

You have three restriction enzymes available, *ApaI*, *EcoRI*, and *SmaI*. Tap the one of these three enzymes that will allow you to distinguish between the four possible plasmids. Please note, that DNA fragments with a size smaller than 100 bp will result in very faint bands and that their size cannot be determined accurately.

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DIGESTION SOLUTIONS (4 POINTS)

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Node Id: 6c455f6de8153f592c46f58b

Restriction enzyme digestion of the two plasmids

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Protocol

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1. Label a strip of 0.2 mL PCR-tubes with S1, S2, S3 and S4 respectively.

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2. Design your digestion experiment to identify the two plasmids using one of the enzymes for all digestions. Include a control with uncut plasmid DNA. The total volume for each reaction should be 10 μL . A typical restriction enzyme digestion includes the following: 2 μL plasmid DNA, 1 μL of restriction enzyme, 1 μL Buffer and water to a final volume of 10 μL .

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Indicate with integers, what you add to each of the four tubes (all amounts are in μL). Write '0' if you do not want to add this specific ingredient

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Q. 4

REACTION CONDITIONS FOR DIGESTION (2 POINTS)

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Node Id: d753f5302f0bf03853b72487

Protocol continued

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3. Mix the necessary components for the restriction digestion reactions as you specified in Question 2.

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4. Set the volume on the pipette to 5 microliter and mix thoroughly by pipetting up and down 5-10 times in each tube. Avoid creating bubbles by keeping the pipette tip below the surface of the liquid at all times.

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5. Place the strip of 0.2 ml tubes in the small Zip lock bag. Close the bag and stick on label with your ID number to the outside of the bag

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Indicate here, the reaction conditions you use for the digestion.

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2. 37°C

Protocol continued

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6. Raise your pink card and an Official will bring your samples to an incubator of your choice.

Complete digestion of the plasmid DNA will require 15 minutes of incubation

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7. After 20 min raise your pink card and ask an Official to return your samples. Your samples will **NOT** be returned unless you ask for it.

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8. Check that you have received your own samples.

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9. Add 3 μ L of loading buffer to each of the samples

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10. Mix by pipetting up and down 10 times (avoid making bubbles)

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Use the Agarose gel electrophoresis' protocol, and analyze the restriction enzyme digestions by gel electrophoresis.

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AGAROSE GEL ELECTROPHORESIS

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WARNING: ELECTRICAL HAZARD – Do NOT move the RunOne electrophoresis apparatus once it has been started! Do NOT insert objects through the lamella in the lid of the apparatus when it has been started!

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WARNING: Potentially HAZARDOUS COMPOUNDS. ALWAYS wear gloves when handling the agarose gel as it is prestained with a DNA binding dye.

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Suggested workflow:

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Your samples from Part 1 and/or 2 should be analyzed by agarose gel electrophoresis.

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Electrophoresis Protocol

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1. Put on rubber gloves.
 2. Take the lid off the RunOne electrophoresis chamber by lifting it straight up.
 3. Prepare the agarose gel for electrophoresis by unpacking it from the white bag. Remove the transparent lid and carefully lift tray with the gel out of the transparent casing. Handle the gel carefully as it easily breaks. Keep the white bag for later use.
1. 戴上塑膠手套
2. 將電泳槽系統的上蓋取下(垂直向上)
3. 打開裝膠的白色塑膠袋，移開蓋子，小心地將膠片盤移出鑄膠盒，不要弄破膠片，保留白色塑膠袋，等一下有用
-
4. Examine the gel for visible damage such as cracks, missing wells or bubbles trapped inside the gel. If you identify such problems report it immediately to the Officials by raising your pink card.
 5. Place the gel as shown in Fig. 1.2.
4. 檢查膠片是否有損傷(裂開、氣泡)，如果發現有缺損，立刻舉粉紅卡向試務人員報告
5. 如Fig. 1.2. 放置膠片


Figure 1.2: RunOne electrophoresis chamber (left) and possible orientations of the agarose gel.
Figure 1.2: 電泳槽系統(左)和電泳膠片放置方向

Electrophoresis protocol continued

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7. The gel should be completely submerged in the running buffer. **If this is NOT the case** please raise your pink card and an Official will help you.
7. 膠片必須完全浸在電泳緩衝液面以下，如果不是，立刻舉粉紅卡，試務人員會協助
8. Load the agarose gel as described in Table 1.2 (samples from Part 1) and Table 1.3 (samples from Part 2). Load 3 μL of 1 Kb ladder, 13 μL sample from Part 1 and 15 μL from Part 2.
8. 依照Table 1. 2 (第一部分樣品)和Table 1. 3的描述，注入3 μL 的1Kb尺標DNA，13 μL 的樣品 (第一部分樣品)，15 μL 的樣品 (第二部分樣品)

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Electrophoresis protocol continued

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9. Place the lid on the RunOne apparatus (the tab slide down into the power supply box).

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10. Check that the voltage is set to 100 V otherwise use to "Voltage select" button.

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11. Switch on the power by pushing the "Run/Stop" button and run the gel for 15 to 20 minutes.

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Note that you will not be able to see the DNA bands migrating through the gel in real time, but that the blue indicator dye in the loading buffer migrates at the same speed as a 300 bp DNA fragment.

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12. After 20 minutes switch off the power by pushing the "Run/Stop" button

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13. Hold up your pink card and an Official will bring you the equipment required for viewing the DNA in the agarose gel and for documenting this.

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Document the DNA band pattern in the gel by taking a photo of it using your tablet and the equipment the IBO Official brings you. The Official will

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14. Put on a pair of gloves.

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15. Transfer the agarose gel in its tray onto the light-table.

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16. Place the black photo hood on top of the light-table.

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17. Place the tablet on top of the photo chamber in a way that the lens is pointing down into the chamber.

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Take three favorable pictures of your DNA gel. These images will be used to judge the success of your PCR and restriction enzyme digestion reactions, loading of a ladder and the validity of the analysis.

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Once you have documented your gel, please transfer it to the white bag, in which it was delivered

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Q. 6

DIGESTION OF PLASMID DNA BY RESTRICTION ENZYME

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Node Id: 631ab61f0da876d8b070db01

Indicate here for each sample, if the plasmid DNA has been digested by the added restriction enzyme, or not.

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
You will now be provided with a standardized gel picture showing results obtained for the restriction enzyme you chose (Question 3). Use this picture to answer Questions 7 and 8. **NOTE YOU WILL ONLY HAVE ACCESS TO THIS PICTURE ONCE YOU HAVE DOCUMENTED YOUR GEL.**

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1 kb plus 尺標DNA

Q. 7

DNA FRAGMENT SIZE(S) OF YOUR RESTRICTION
ENZYME DIGESTION

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Node Id: 723065e993cf22c5ad91aad4

Indicate here, with a tag if a DNA fragment of a given size was present in
your restriction enzyme digestion

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Q. 8

IDENTIFICATION OF PLASMIDS IN TUBES

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Node Id: f182a828e66bfc5be092fde

Indicate here for each tube which plasmids are present, and which plasmids are absent. Indicate with a '+' in the last column if it was impossible to determine the presence of plasmids for a tube.

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PLASMID
PZ

PLASMID
PZ:ALB

IMPOSSIBLE
TO
DETERMINE
BASED ON
THE CHOSEN
ENZYME
依選擇的酶

Plasmid tube 1	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>
Plasmid tube 2	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>

2. PCR BASED GENOTYPING OF YEAST MUTANTS (37 POINTS)

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A mutant yeast strain was isolated. This yeast strain, requires tyrosine and phenylalanine in the growth medium to be able to propagate. The ability of the mutant to grow on a medium without tryptophan has not been tested. A back cross of the mutant strain has shown that only a single gene is affected.

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Your task is now to determine the genetic basis for the observed auxotrophy. You have five primer pairs available (Table 2.1), each capable of detecting dysfunctional (mutant) alleles of the genes encoding five key enzymes in the biosynthetic pathway for formation of aromatic amino acids (Fig. 2.1).

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Table 2.1. The available primer pairs, the genes they amplify and the expected DNA amplicon in the wild type and mutant alleles.

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
Figure 2.1: Biosynthetic pathway for the formation of aromatic amino acids (from upper left: Tryptophan, Tyrosine, Phenylalanine, Prephenate, Chorismate and Shikimate) and genes encoding the individual enzymes in the pathways.

Second Primer Pair 第二個引子對

Protocol

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1. Make a Master-mix for the PCR reactions in a 1.5 ml tube by mixing: 196 μ l Water + 70 μ l DNA polymerase buffer (x5) + 35 μ l dNTP's + 7 μ l DNA polymerase.

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2. Mix by pipetting up and down 5-10 times with the pipette set at 100 μ l.

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3. Find the strip of 0.2 ml PCR tubes and glue on a small label with your ID number.

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4. Label the tubes 1 - 6 (**DO NOT** break the strip of tubes into individual tubes – if you do so, your samples will not be accepted for PCR).

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5. Pipette 44 μ l of the Master-mix to each of the six PCR tubes.

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6. To tubes 1 to 3, add 5 μ l of the first Primer pair you have chosen.

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7. To tubes 4 to 6, add 5 μ l of the second Primer pair you have chosen.

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11. Mix the individual reaction by pipetting up and down 5-10 times with the pipette set at 45 μ l.

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12. Close the lids.

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Table 2.2. Three PCR programs (1-3), stored at 10°C.

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Program 1 PCR反應程式1	Program 2 PCR反應程式2	Program 3 PCR反應程式3
98°C for 1 min	98°C for 1 min	98°C for 1 min
25 times: 98°C for 15 sec 55°C for 5 sec 72°C for 20 sec	25 times: 72°C for 20 sec 55°C for 5 sec 98°C for 15 sec	25 times: 72°C for 10 sec 98°C for 15 sec 55°C for 5 sec
72°C for 1 min	72°C for 1 min	72°C for 1 min
Store at 10°C	Store at 10°C	Store at 10°C

Q. 10 CHOICE OF PCR PROGRAM (2 POINTS)

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Node Id: 33418d50512f939d3aeb6813

Indicate here, which of the three PCR programs (Table 2.2) should be used for successful amplification if the primers have a melting temperature (T_m) of 60 °C and the used DNA polymerase can synthesize DNA at a speed of 25 bp/sec at 72 °C?

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3. PCR PROGRAM 3 PCR反應程式3

IMPORTANT

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1. Assure yourself that Questions 9 and 10 are answered as you intended. Once you are happy with the answers, lock them by clicking "LOCK ANSWERS".

Note: after locking, you will not be able to change your answers anymore!

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2. Once your answers are locked, raise your pink card and an official will transfer your samples to the thermal cyclers. The PCR program takes 20 minutes to complete. (In the meantime go to Parts 1 or 3 according to the time table in the introduction).

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3. After 20 minutes raise your pink card and an official will return your samples. Please note that the sample will NOT be returned unless you ask for it.

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4. Upon receiving your samples please check that they are **your** samples.

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5. Prepare the PCR samples for agarose gel electrophoresis by adding 10 µL of loading buffer to each of the samples and mix by pipetting up and down 5 times (avoid making bubbles).

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6. Load your samples as described in the Agarose gel electrophoresis section.

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7. Documentate the gel by taking a photo (the photo of this gel will be inspected and a maximum of 22 points will be awarded)

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You will now be provided with a standardized gel picture based on the

primer pairs you previously chose to use for the PCR analysis. use this picture to answer Questions 11 and 12. **NOTE: YOU WILL GET THIS PICTURE FROM THE ASSISTANT RETURNING YOUR GEL.**

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Q. 11 SIZE OF DNA BANDS (4 POINTS)
DNA 條帶(bands)大小

Node Id: 2261c610758917fbc24905ed

Indicate here if a given size of the obtained DNA (bp) products was present '+'. Report your results column-by-column.

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Q. 12 DYSFUNCTIONAL ENZYME STEPS (4 POINTS)
失效酵素作用步驟(4分)

Node Id: ff6dc825d5bde0184dfec782

Based on your results, indicate for each of the five possible enzymatic steps shown in Fig. 2.1 if they are functional or dysfunctional, or whether their status can not be determined with the used primers.

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3. AMINO ACID AUXOTROPHY OF MUTANT YEAST (15 POINTS)

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For many years, the haploid fungus *Saccharomyces cerevisiae* (baker's yeast) has been used as a model organism for elucidating the central metabolism in eukaryotes. Yeast is normally capable of synthesizing all twenty amino acids needed for protein synthesis. Following UV-mutagenesis various amino acid auxotrophic mutants have been isolated. 多年來，單倍體麵包酵母菌被用來做為解答真核生物主要代謝途徑的模式生物。正常酵母菌能自行合成所有20種胺基酸。利用紫外光(UV)突變誘導後，分離出各種不同胺基酸合成缺失突變株。

The booklet with eight photos depicts haploid *S. cerevisiae* strains 1–5 grown on eight different media (A to H). 照片冊中的8張照片呈現了5個單倍體酵母菌株(1~5)在8種培養基(A~H)上的生長情形：

Position 1 = no inoculum, Position 2 = Strain 1; Position 3 = Strain 2; Position 4 = Strain 3; Position 5 = Strain 4; Position 6 = Strain 5. The media composition for plate A to H is given in Table 3.1

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Q. 13

GROWTH/LACK OF GROWTH OF YEAST STRAINS (10 POINTS)

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Node Id: bbd24126db8bb263c5dc6452

Record here growth or lack of growth for the various strains, use: "+" for growth and "-" for no growth.

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Node Id: 5d4effda2f48e112ec5b7354

Based on the recorded growth patterns, deduce which enzymatic step(s) that is/are most likely dysfunctional, if any, in the five strains (1–5). For each mutant write the digit (1–31) for the dysfunctional enzymatic step (see Biochemical chart on the wall of the workspace) or 0 if no step is dysfunctional.

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	STRAIN 1 菌株1 STRAIN 1	STRAIN 2 菌株2 STRAIN 2	STRAIN 3 菌株3 STRAIN 3	STRAIN 4 菌株4 STRAIN 4	STRAIN 5 菌株5 STRAIN 5

ENZYME STEP NO. in biochemical chart on the wall 酵素作用步驟					
--	--	--	--	--	--


Figure 3.1: Simplified model for amino acid biosynthesis in yeast.

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