Country:

Student Code:

27th International Biology Olympiad

17th-23rd July, 2016 Hanoi, Vietnam



Practical Exam 3 **BIOCHEMISTRY AND MICROBIOLOGY Answer Key**

Total points: 100


EXPERIMENT 1 (40 points)

Q.1.1. (6 points)

	Tube 1	Tube 2	Tube 3	Tube 4	Tube 5	Tube 6
Sample	Protein Marker	No_IPTG	IPTG	Pellet	Super	Puri-P
Stock (µl)	5	10.4	17.9	14.9	14.9	15.0
H ₂ O (μl)	5	14.8	11.0 - 11.1	12.5	12.5	12.5
2X SDS-PAGE loading buffer (μl)	10	14.8	11.0 - 11.1	12.6	12.6	12.5
Total volume (μl)	20	40	40	40	40	40

^{0.5} point for each correct answer for Stock.

Q.1.2. (10 points)




5 points for SDS-PAGE picture with bands of sample.

Add I point for clearly visible bands in each well except Marker.

^{0.4} point for each correct answer for H_2O and loading buffer.

Q.1.3. (4 points)



2 points for at least 5 correct spots on the graph.

2 points for the linear curve.

Q.1.4. (4 points)

Protein	Size (kDa)
H	20 ± 2
В	35 ± 2

Scoring bases on the student's answer in Q.1.3.

2 points for each correct answer.

Q.1.5. (4 points)

Statement	True	False
A. H protein is over expressed in LB media with IPTG.	✓	
B. B protein is completely soluble in the nickel binding buffer.		✓

C. H and B protein interact with each other.D. Majority of recombinant proteins bound to the nickel		✓	
D. Majority of recombinant proteins bound column.	to the nickel	✓	

⁴ points for 4 correct answers.

2 points for 3 correct answers.

1 points for 2 correct answers.

0 points for 0-1 correct answers.

Q.1.6. (4 points)

	Statement	True	False
A.	SalI and BamHI can be used to insert b gene into MCS-1	✓	
В.	Gene <i>h</i> and <i>b</i> should be cloned in the same orientation to be expressed simultaneously.	√	
C.	Gene h and b should be in the same reading frame to be expressed simultaneously.		/
D.	To maintain the plasmid, ampicillins should be added to the culture medium.	√	

⁴ points for 4 correct answers.

0 points for 0-1 correct answers.

Q.1.7. (4 points)

Sample	Size (kDa)
The peak observed in sample 1	50 ± 1,0
The peak observed in sample 2	100 ± 20

Q.1.8. (4 points)

Statement	True	False
A. H protein exists as monomer.		✓
B. H and B probably exist as heterodimer.		✓

² points for 3 correct answers.

¹ points for 2 correct answers.

C. H protein helps to stabilize B protein D. In native gel-filtration column analysis, retention time of a		
D. In native gel-filtration column analysis, retention time of a		
protein is proportional to its molecular weight.		· · · · · · · · · · · · · · · · · · ·

Scoring bases on the student's answer in Q.1.7.

⁴ points for 4 correct answers.

² points for 3 correct answers.

I points for 2 correct answers.

⁰ points for 0-1 correct answers.

EXPERIMENT 2. (30 points)

Q.2.1. (4 points)

Ascorbic acid sample				Coffee ex	tract sample			
Diluted	Volume	Volume	Concentration	Diluted	Volume	Volume	Concentration	
solution	(µl) of	of H ₂ O	(mg/ml)	Solution	(µl) of	of H ₂ O	(mg/ml)	
	solution	(µl)			solution	(µl)		
	used				used			
	for				for			
	dilution				dilution			
AA1	40	160	0.2	CC1	200	0	5	
AA2	100	100	0.1	CC2	100	100	2.5	
AA3	100	100	0.05	CC3	100	100	1.25	
AA4	100	100	0.025	CC4	100	100	0.625	

^{0.25} point for each correct answer.


Q.2.2. (5 points)

Ascorbic acid				Coffee extract				
Solution	tion Log ₁₀ of Mean absorbanc (mg/mL)		SC%	Solution	Log ₁₀ of concentration (mg/mL)	Mean absorbance	SC%	
Control		0.90 ± 0.01						
AA1	- 0.70	0.076 ± 0.007	91.48 ± 9.1	CC1	0.70	0.225 ± 0.022	75.01± 7.5	
AA2	- 1	0.249 ± 0.024	72.34 ± 7.2	CC2	0.40	0.334 ± 0.033	62.88 ± 6.2	
AA3	- 1.30	0.511 ± 0.051	43.15 ± 4.3	CC3	0.10	0.602 ± 0.060	33.06 ± 3.3	
AA4	- 1.60	0.716 ± 0.071	28.10 ± 2.8	CC4	- 0.20	0.744 ± 0.074	17.26 ± 1.7	

^{0.2} point for each correct answer.

Scoring values of SC% based on the value of mean absorbance.

Q.2.3. (5 points)



Log₁₀ concentation (mg/mL)

3 points for 4 correct spots on the graph.

2 points for the linear curve.

Q.2.4. (5 points)

For example: SC₅₀ of ascorbic acid


$$Log(SC_{50}) = -1.23$$

$$SC_{50} = 10^{-1.23} = 0.058$$

Ascorbic acid	Coffee extract
0.06 ± 0.006	2.01 ± 0.20

Scoring values of SC_{50} based on the graph.

Q.2.5. (3 points)



Q.2.6. (4 points)

Statement	True	False
A. Antioxidant activities of the diluted coffee extracts are negligible.		✓
B. To obtain more accurate determination of antioxidant activities, another experiment with higher concentration samples needs to be carried out.		✓
C. Activity of antioxidant enzymes in coffee extract resulted in the above result.		√
D. If NADH is added in the wells, there will be no change in the assumed absorbance	√	

⁴ points for 4 correct answers.

I points for 2 correct answers.

0 points for 0-1 correct answers.

Q.2.7. (4 points)

Statement	True	False
A. The student has used DPPH 0.1 mM for Protocol A*.		✓
B. The student has loaded 10 µl of sample in each well for Protocol A*.	√	
C. After adding DPPH, the student has incubated the 96-well for a shorter time than in Protocol A.	√	
D. The student has used better solvent for antioxidants.		V

⁴ points for 4 correct answers.

² points for 3 correct answers.

² points for 3 correct answers.

I points for 2 correct answers.

⁰ points for 0-1 correct answers.

EXPERIMENT 3. (30 points)

Q.3.1. (3 points) $2 \pm 0.1 \text{ h}$

Q.3.2. (3 points) $4.2 \times 2 \times 10^8 \times 2^3 = 67.2 \times 10^8$ cells

Q.3.3. (2 points)

Cultivation time (h)	Sample	Dilution factor	Sample volume (mL)	Deionized water volume (mL)	Total volume (mL)
0	A0	2	15	15	30
6	A2	5	6	24	30
9	A3	10	3	27	30
15	A5	20	1.5	28.5	30

0.25 point for each correct answer.

Q.3.4. (9 points)

Cultivation time (h)	Sample	Volume of 0.1 M NaOH titrated	Volume of 0.1 M NaOH titrated	Volume of 0.1 M NaOH titrated
		(First time, mL)	(Second time, mL)	(Average, mL)
0	A0	3.8 ± 0.1	3.8 ± 0.1	3.8 ± 0.1
6	A2	2.2 ± 0.05	2.2 ± 0.05	2.2 ± 0.05
9	A3	2.4 ± 0.05	2.4 ± 0.05	2.4 ± 0.05
15	A5	2.7 ± 0.05	2.7 ± 0.05	2.7 ± 0.05

0.75 point for each correct answer.

Q.3.5. (10 points)

Cultivation time (h)	Sample	Volume of 0.1M NaOH (mL) need to be used to titrate 30 mL of stock sample	Lactic acid concentration (g/L)
0	A0	7.6 ± 0.2	0
6	A2	11 ± 0.25	(A2-A0)*9/30 = 1.02
9	A3	24 ± 0.5	(A3-A0)*9/30 = 4.92
15	A5	54 ± 1	(A5-A0)*9/30 = 13.92

1 point for each correct answer for Volume of 0.1 M NaOH (scoring based on the result of Q.3.4).

2 points for each correct answer for lactic acid concentration.

Q.3.8 (3 points) $0.33 \times 2 \times 10^8 + 6 \times 2 \times 10^8 = 12.66 \times 10^8$ tế bào.